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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Low elastic orthosis with rigid and thermoformable pads to stabilize the lumbar rachis and relieve pressure

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF. 1078 · 1078N					
Size	XS	S	M	L	XL	XXL
Circumf. of waist cm	50/60	60/75	75/90	90/105	105/120	120/135
Front height cm	20	20	20	20	20	20
Back height cm	26	26	26	26	26	26
Colour	1078 = beige · 1078N = black					

OPTION

9ST012/L – optional rigid stays 22 cm

MAINTENANCE

 Do not bleach  No chemical cleaning
 Do not iron  Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Dispose of this product and its parts responsibly

INDICATIONS

- Backache
- Minor injuries to the lumbosacral spine and paravertebral lumbar contractures

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Cotton and elastane fabric
- Front Velcro® fastening system
- Back reinforced by crossed straps
- 2 front and back thermoformable pads, on medical prescription
- Optional rigid stays, mouldable on medical prescription; to insert in the 4 linear back pockets
- Adapts to the pathology of the patient

PUTTING ON THE APPLIANCE

- Open the corset and position it so that the loop (A) is at the top.
- Place a thumb in the loop (A) while donning the corset so that it easier to close and fasten securely with the central Velcro® fastening (B). Ensure that the fastening is properly centred.

NOTE:

- To apply the rear rigid pad (heat-mouldable):
- apply the adhesive padding to the plastic part (concave side-without any writing, fig.C)
 - apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig.D)
 - attach the pad with Velcro® to the internal dorsal part of the orthosis

- To apply the front rigid pad (heat-mouldable):
- apply the adhesive padding to the plastic part (concave side-without any writing)
 - attach the pad with Velcro® to the front closure of the orthosis

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